



Anti-E2 Tag Mouse Monoclonal Antibody(3E7)

Catalog No.	EM33016-01	EM33016-02
Size	50 μ l	100 μ l

E2-tag peptide GVSSTSSDFRDR conjugated to KLH. E2- Tag antibody can recognize C-terminal, internal, and N-terminal E2-tagged proteins.

Molecular Weight: N/A


Application : WB 1:5000 Optimal dilutions should be determined by the end user.

Specificity:The E2 tag antibody can recognize E2 tag fusion proteins.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



2 μ g E2 fusion protein+ Primary antibody dilution at 1、 1:5,000 2、 1:10,000