



## Anti-HSV Tag Mouse Monoclonal Antibody(8C5)

Catalog No.	EM33017-01	EM33017-02
Size	50 $\mu$ l	100 $\mu$ l

HSV-tag peptide QPELAPEDPED conjugated to KLH. HSV- Tag antibody can recognize C-terminal, internal, and N-terminal HSV-tagged proteins.

**Molecular Weight:** N/A

**Application :** WB 1:5000 Optimal dilutions should be determined by the end user.

**Specificity:**The HSV tag antibody can recognize HSV tag fusion proteins.

**Form:** Liquid ,1mg/ml

**Storage buffer:** PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

**Storage:** Store at -20°C. Do not aliquot the antibody.

