



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-Avi Tag Mouse Monoclonal Antibody(5D2)

Catalog No.	EM33036-01	EM33036-02
Size	50 μ l	100 μ l

The Avi tag is a biotin-acceptor peptide¹, GLNDIFEAQKIEWHE. The 15-residue peptide served as a substrate mimic for biotin ligase (BirA), which usually recognizes the much larger protein domain.

Molecular Weight: N/A


Application : WB 1:3,000-5,000 Optimal dilutions should be determined by the end user.

Specificity:The Avi-Tag Mouse Monoclonal Antibody can detects Avi-Tag recombinant protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of Avi-Recombinant protein with AVI mAb diluted at
1) 1:5,000 2) 1:10,000

EMAR

EMAR