



Beijing Emarbio Science & Technology Co.,LTD

## Anti-RFP Mouse Polyclonal Antibody

Catalog No.	EM34020-01	EM34020-02
Size	50μl	100μl

The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.

**Molecular Weight:** N/A


**Application :** WB 1:5,000 IHC 1:200 IP 1:100 Optimal dilutions should be determined by the end user.

**Specificity:** The RFP tag antibody can recognize RFP tag fusion proteins.

**Form:** Liquid ,1mg/ml

**Storage buffer:** PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

**Storage:** Store at -20°C. Do not aliquot the antibody.



2ug RFP fusion protein+ Primary antibody dilution at  
1、 1:5,000    2、 1:10,000