



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-ECFP Tag Mouse Polyclonal Antibody

Catalog No.	EM34024-01	EM34024-02
Size	50 μ l	100 μ l

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Molecular Weight: N/A


Application : WB 1:5,000 IHC 1:200 IF 1:200 Optimal dilutions should be determined by the end user.

Specificity:The ECFP Antibody can detects ECFP & CFP tag protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of ECFP Recombinant protein with ECFP PAB diluted at
1) 1:5,000 2) 1:10,000

EMAR

EMAR