



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-mStrawberry Tag Mouse Polyclonal Antibody

Catalog No.	EM34027-01	EM34027-02
Size	50µl	100µl

mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien' s lab by directed mutagenesis of mRFP.

Molecular Weight: N/A

Application : WB 1:5,000 IHC 1:200 IF 1:200 IP 1:200

Optimal dilutions should be determined by the end user.

Specificity:The mStrawberry Mouse Polyclonal Antibody can detects mStrawberry & mStrawberry Tag protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of mStawberry recombinant protein with mStrawberry PAB diluted at 1:5,000.

EMAR

EMAR