



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-ECFP Tag Rabbit Polyclonal Antibody

Catalog No.	EM34029-01	EM34029-02
Size	50 μ l	100 μ l

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Molecular Weight: N/A

Application : WB 1:2,000-5,000 IP 1:200

Optimal dilutions should be determined by the end user.

Specificity: The ECFP Antibody can detect ECFP & CFP tag protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of ECFP recombinant protein with ECFP PAB diluted at 1:5,000.

EMAR

EMAR