



EMAR Beijing Emarbio Science & Technology Co.,LTD

## Anti-mStrawberry Tag Rabbit Polyclonal Antibody

Catalog No.	EM34032-01	EM34032-02
Size	50µl	100µl

mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien' s lab by directed mutagenesis of mRFP.

**Molecular Weight:** N/A

**Application :** WB 1:2,000-5,000 IHC 1:200 IF1:200 IP 1:200

Optimal dilutions should be determined by the end user.

**Specificity:**The mStrawberry Mouse Polyclonal Antibody can detects mStrawberry & mStrawberry Tag protein.

**Form:** Liquid ,1mg/ml

**Storage buffer:** PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

**Storage:** Store at -20°C. Do not aliquot the antibody.



Western blot analysis of mStawberry recombinant protein with mStrawberry PAB diluted at 1:5,000.

EMAR

EMAR